

Laerdal Airway Management Trainer

The skills trainer for adult airway management training

Realistic practice is the key to developing proficiency in airway management skills. The Laerdal Airway Management Trainer's lifelike upper torso and head simulates real-world complications when practicing a variety of intubation, ventilation, and suction techniques.



Product features:

- Practicing of oral and nasal intubation
- Practicing use of LMA (Laryngeal Mask Airway) and Combitube
- Correct tube placement can be checked by practical inflation test
- Realistic anatomical features allow demonstration of Sellick Maneuver and laryngospasm
- Bag-Valve-Mask ventilation can be practiced
- Stomach inflation and vomiting situation can be simulated
- Provides visual inspection of lung expansion
- Provides auscultation of breath sounds
- Airway demonstration model is standard with each trainer

25 00 00 Laerdal Airway Management Trainer

Laerdal Airway Management Trainer Includes:

Manikin on Sturdy Board, Airway Demonstration Model, Cleaning Kit, Airway Lubricant, Carry Case, and Directions for Use.

Product benefits:

- Realistic representation of human anatomy, tissue, and skin
- Allows students to undertake training that is directly transferable to the clinical setting
- Practical training in clearing an obstructed airway and suctioning of liquid foreign matter
- Minimal maintenance and robust design delivers cost-effective training
- Base plate mount allows stable practicing conditions
- A hard carry case provides easy transportation and safe storage

Consumables



25 28 00
Concentrated
Simulated Vomit



25 20 90
Airway Lubricant