

Laerdal ALS Baby

The realistic skill trainer for practicing training in infant emergencies

The Laerdal ALS baby is a portable skill trainer for realistic infant resuscitation training. Laerdal ALS Baby manikin is a three-month-old infant with an anatomy that offers exceptional realism for individual training. The Laerdal ALS baby offers the opportunity to practice advanced resuscitation skills, including airway management, professional rescuer CPR, vascular access, and 4-lead ECG monitoring.



ALS Baby is compatible with the HeartSim 4000 Rhythm Simulator. Must be purchased separately, Computer not included (call for further details)

Product benefits:

- Educational effectiveness - through practicing basic and advanced CPR skills
- Multifunction use - facilitates training of multiple skills needed for infant ACLS training
- Cost efficient - durable and dependable for long-term use and cost efficiency
- Anatomically realistic enables lifelike infant resuscitation training to be performed
- Portable infant ACLS skill trainer for training in all environments

Product features:

- Realistic airway anatomy with tongue, oropharynx, epiglottis, larynx, vocal cords and trachea
- practicing of bag-valve-mask ventilation, oral and nasal intubation and use of LMA (Laryngeal Mask Airway)
- CPR can be performed
- 3-lead, 4 connectors ECG monitoring
- Available with HeartSim 200 Rhythm Simulator
- Intraosseous needle insertion with aspiration of bone marrow
- Sellick Maneuver can be performed



HeartSim 200 Cardiac Rhythm Simulator can be purchased separately.

08 00 31 Laerdal ALS Baby 200 (Complete)

08 00 30 Laerdal ALS Baby (manikin only)

26 00 10 HeartSim 200 Cardiac Rhythm Simulator

Laerdal ALS Baby 200 (Complete) Includes:

ALS Baby Manikin, Molded Plastic Carry Case, 5 IO Leg Replacement Pads, Pulse Bulb, Lubricant Jelly, HeartSim 200 and Directions for Use.

Laerdal ALS Baby Includes:

ALS Baby Manikin, Molded Plastic Carry Case, 5 IO Leg Replacement Pads, Pulse Bulb, Lubricant Jelly and Directions for Use.

Consumables



25 20 90
Airway Lubricant



08 23 05
Leg Replacement
Pads (pkg. 5)



08 20 20
Lungs/Stomach
(pkg. 6)