


*to protect both the patient and the rescuer.*



## For all your CPR Barrier Device needs...

*It's up to every family member, co-worker, security guard, school nurse, friend and stranger to step in and take the appropriate action to help save a life before the ambulance gets there!*

*And it's Laerdal's mission to help save those lives by providing the high quality products that continue to improve the quality of CPR. Ever since 1969 when Laerdal introduced the very first folding Pocket Mask, we have created, supplied, and improved CPR solutions that protect and help both the patient and the rescuer.*

*We have everything you need to help give patients the best possible CPR. That's why all you need is Laerdal for all your CPR Barrier Device needs.*

[www.vascomedical.gr](http://www.vascomedical.gr)



**Laerdal**<sup>®</sup>  
helping save lives

## Laerdal Pocket Mask

*The Laerdal Pocket Mask is a world-leading CPR mask!*

The Laerdal Pocket Mask incorporates a one-way valve and 3M™ Filtrete™ fi lter allowing for effective ventilation and outstanding user protection. While there may be look-a-likes on the market, there is only one real Laerdal Pocket Mask.!



## Laerdal Paediatric Pocket Mask

*The most widely used CPR mask in the world is now designed with infants and children in mind.*

Utilizing the same low-resistance one-way valve as the traditional Laerdal Pocket Mask in combination with a transparent soft silicone mask, the Laerdal Paediatric Pocket Mask is an excellent companion for the professional rescuer and fi rst responder who have contact with infants and children.



## Laerdal Face Shield

*The Laerdal Face Shield with 3M™ Filtrete™ fi lter is small in size, but BIG in protection for user and victim!*

The Laerdal Face Shield for single use is available in convenient foil packages or sealed key rings.



460000 Laerdal Face Shield CPR Barrier

## Laerdal Face Shield Key Ring

*The Laerdal Face Shield with 3M™ Filtrete™ fi lter is small in size, but BIG in protection for user and victim!*



460008 Laerdal Face Shield CPR Barrier Key Ring



**Vasco Medical**

36A Taxiarchon str. 17455 Alimos

Tel +30 210 9834733

Fax +30 210 9847016