



**Laerdal®**  
helping save lives



## Quality CPR Training with Real-Time Feedback

### Resusci Anne® Wireless SkillReporter™

The new Resusci Anne® Wireless Skill Reporter™ is the latest innovation in the Laerdal portfolio. Pairing the RA SkillReporter manikin with an easy-to-use software program, the Resusci Anne Wireless SkillReporter provides real-time feedback throughout skills practice and a debriefing screen to evaluate a student's performance – building real CPR competencies during training.

The Resusci Anne family of manikins set the standard for quality CPR manikins, now with the new Resusci Anne Wireless SkillReporter, Laerdal is setting a new standard for CPR skills practice.



laerdal.com



# Resusci Anne® Wireless SkillReporter™

## FEATURES

- Simple-to-use software
- Enhanced user interface
- Multiple-manikin class overview (up to 6 manikins)
- Real-time performance feedback
- Comprehensive debriefing tools
- Software compatible with current Resusci Anne® SkillReporter Manikin

## SPECIFICATIONS

### Minimum Software Requirements

- Windows XP, Windows Vista or Windows 7

### Minimum Hardware Requirements

- Core 2 Duo or better
- 1 GB RAM (2 GB recommended)
- 200 MB free hard disk space
- 1024 x 768, 1280 x 800, 1280 x 1024 or better
- 16-bit color resolution or better
- Optical drive required for installation
- Internet connection for automatic license activation
- Bluetooth® required for operation

## ORDER INFORMATION

**318000 Resusci Anne® Wireless SkillReporter™ Complete\***  
Includes: Software, Wireless Kit and Bluetooth® Adapter  
\* Manikin and Laptop PC not included

**31005601 Resusci Anne® SkillReporter™ Manikin**

## ACCESSORIES

**318001 Wireless SkillReporter™ USB Adapter**  
Includes: USB Bluetooth Antenna

**318002 Resusci Anne® Wireless Kit**  
Includes: Bluetooth Serial Adapter and Manikin Cable Adapter

**318003 Resusci Anne® Wireless SkillReporter™ Software**  
Includes: Software CD only

Performance  
Debrief



Multi-Manikin  
Class Overview



Real-Time  
Feedback



See Resusci Anne  
Wireless SkillReporter  
in action



877-LAERDAL (523-7325)  
laerdal.com



**Laerdal®**  
helping save lives