



Laerdal[®]
helping save lives



More Efficient CPR Training

Resusci Anne[®] QCPR

The Resusci Anne[®] QCPR manikin and technology enables healthcare providers to obtain competence to deliver quality CPR in accordance with the latest Guidelines. QCPR technology is utilized to help both instructors and trainees more effectively monitor and review CPR performance, increasing the efficiency and value of training time.



Resusci Anne® QCPR

In addition to providing high-quality CPR training, the Resusci Anne® QCPR manikin allows for those skills and others to be integrated into team training. The available add-on options and accessories make Resusci Anne QCPR more than just a high-quality CPR training manikin. The manikin can be upgraded to train teams on defibrillator use, airway management, IV administration, and rescue/trauma situations.

Key Training Features

CPR

- Quality CPR feedback (with feedback device)
- Mouth-to-mouth ventilations
- Ventilations with face shield, pocket mask, and BVM

Airway

- Realistic chest rise/fall with ventilations
- Sellick maneuver
- Positive pressure ventilation
- Airway occlusion
- Supraglottic airway (with optional Airway Head)

Pulses

- Manual Carotid pulses

Defibrillation

- Train using a live defibrillator (with ShockLink™)
- Detect Pad Placement
(with optional AED Skin, plus Laerdal AED Trainer or ShockLink)

General

- Interchangeable pupils (with optional Airway Head)
- IV initiation and fluid administration (with optional IV Arm)
- Rescue/Trauma Simulation
(with optional Trauma Limbs or Rescue Limbs)
- Latex free
- Wired and wireless feedback device options
- Bluetooth Low Energy
- Rechargeable battery (up to 34 hours operation time per charge)

QCPR Feedback

- Compression rate and depth
- Complete release
- Correct hand position (not available on apps)
- Appropriate ventilation volume
- Amount of compressions and ventilations
- Frequency and length of interruptions (not available on SkillGuide/apps)

Feedback Device Options



206-30001
SimPad® PLUS with SkillReporter™
206-50050
SkillReporter Software License



170-30050
SkillGuide™



202-56050
SkillReporter (PC)



Feedback Apps (free)

QCPR Learner

QCPR Instructor

Optional Limb Accessories



312025 /
150-10150
IV Arm



090050
Bleeding
Control Legs



312050
First Aid/
Trauma Limbs



312000
Rescue Limbs

Starting Configurations

Resusci Anne® QCPR

| | |
|-----------|--------------------------------|
| 171-00160 | Torso |
| 171-01260 | Full Body |
| 172-00160 | w/ Airway Head – Torso |
| 172-01260 | w/ Airway Head – Full Body |
| 173-00160 | AED - Torso |
| 173-01260 | AED - Full Body |
| 174-00160 | AED w/ Airway Head - Torso |
| 174-01260 | AED w/ Airway Head - Full Body |

Kits

| | |
|-----------|---|
| 173-50750 | AED Skin Upgrade Kit (add AED skin to non-AED skin configurations) |
| 171-20050 | Periodic Maintenance Kit |
| 316000 | Airway Trainer Update Kit (add Airway Head to non-Airway Head configurations) |

Warranty

One-year limited warranty. Refer to Laerdal Global Warranty for terms and conditions.

Information is subject to change. To learn more, visit Laerdal.com